



MAGNOSEEK

Quick Start Guide

Get started in minutes and experience
smarter wellness today!



Index


Unboxing your MagnoSeek_____	Page 2
Connecting Power_____	Page 3
Bootup & Connection_____	Page 4
Wi-Fi Connection_____	Page 5
Quick Tips for Best Use_____	Page 6
Start Scanning_____	Page 7
Scanning with your Smart Phone_____	Page 8-10
Scanning Without your Smart Phone_____	Page 11-12
Signal Interpretation Chart_____	Page 13
Warranty Coverage_____	Page 14
MagnoSeek Support and User Manual QR _____	Back Cover



Unboxing your MagnoSeek

Verify package contents:

A



MagnoSeek Device

B



**MagnoSeek Wristband
& Attached Cable**

C



Power Adapter

D



Power Cord

E



MagnoSeek Case Keys




Connecting Power

Step 1: Connect the cable (D) to the adapter (C).


Step 2: Plug the charging cable (C) into the device (A).

Tip: Keep the MagnoSeek plugged in while in use.

Step 1



Step 2





Boot Up & Connection

Step 3: Ensure the MagnoSeek Wristband and Cable (B) are attached properly to the MagnoSeek device (A).


Step 4: Press the Power Button, located on the left or right bottom of the MagnoSeek Device (A), until the indicator lights up.

Tip: If the wristband is not connected properly, you will see the following error message


No devices found! Program will terminate!

OK

Step 3



Step 4





Wi-Fi Connection

Step 5: When activated, use the touchscreen to connect the device to your preferred local Wi-Fi.

Tip: If the device is not correctly connected to Wi-Fi, you will receive an error code

Bind Card

Card#	13974100075
Mobile	13974100075
Code	

Fail:

Bind



Quick Tips for Best Use

Please refer to the Official User Manual before operating (See QR code on back cover)

- Not for use by children under 7, pregnant or menstruating women, pacemaker users, or those with neurological conditions.
- For clearer readings, avoid taking medication, coffee, tea, or strong stimulants at least 2 hours before your session.
- Remove glasses (1 hr before), metal, energy stones, and magnetic items.
- Keep phones/electronics at least 1 meter away during operation.
- Do not cross arms or legs; stay relaxed.
- Avoid using other frequency-based devices within 24 hrs before your session.
- Use the restroom beforehand for comfort.
- Store in a cool, dry place.

Start Scanning



Step 6



Place the wrist strap on your left wrist, about one finger width below the palm (at the pulse area), to measure frequencies.

Step 7



Using the touchscreen, select the application/icon "Prife"

From this point, you can scan and view your reports by using 2 options:

Option 1: For a smoother and more streamlined experience, use your smartphone. (Start on page #8)



Option 2: Alternatively, you can use the manual options on the MagnoSeek device. (Start on page #11)




Start Scanning with Smart Phone

Step 8a



Scan the provided QR Code to initiate the registration process and enter your personal details.

Step 9a

A screenshot of the Prife International registration form. The form includes fields for Phone, Name, Birth (8-digit birth date, e.g. 20180808), Gender (M), Blood Type (A), Test Type (Detailed), and AI Repair (0). There is a checkbox for "I have read and agree to User Agreement" and a "Start Smart Detection" button at the bottom.

Fill out the registration: Phone Number, Name (keep it short; no need for full name), Date of Birth (format: YYYYMMDD), Gender, and Blood Type (if applicable).


Step 9a Tips:

If blood type is unknown, select A

Select "Detailed" Test Type

AI Repair: Select 0 just to see your report, select 1-20 to initiate that number of repair functions you would like, after your scan is complete. Each repair takes an estimated time of 40-60 minutes to complete.

If your date of birth is incorrect, you will receive an error notification




Step 10a: Tap "Start Smart Detection"

Step 11a: Wait while MagnoSeek measures and analyzes.

Step 12a: Explore Your Results with Smart Phone



1




After the scan is complete, press the exit symbol to leave Smart Scan.

2



Select "Reporting" at the top right corner to generate your overview report.

3



At the Report Settings window, uncheck any sub-test records you don't need, then press "Collect".

4


Your detailed report will generate.

5

After the report is finished generating, select "AI rpt".




Step 12a: Explore Your Results with Smart Phone



A QR Code will appear labeled "View", scan the code with your smart device.

If this is not your first time using the device, another QR code will appear under "History" which shows all previous scans.




After you select "View" your detailed report will appear on your smart device.

For full explanations of detailed scanning, report generation, energy optimization, system functions, and/or report Interpretation, scan the QR Code on the back cover to access the MagneSeek: Official User Manual.



Start Scanning Without Smart Phone

Step 8b



When you see the "Smart Scan" QR pop-up appear, press the X button to exit and proceed with the manual process.

Select the red journal icon in the top left to initiate the registration process and enter your personal details.

Step 9b

new customers	
surname	
name	
middle name	
birthday	1900/01/01
sex	M
Blood	unknown
Kindred	Unknown
address	
Phone	

Fill out the registration: Phone Number, Name (keep it short; no need for full name), Date of Birth (format: YYYYMMDD), Gender, and Blood Type (if applicable).

Step 9b Tips:

If your date of birth is incorrect, you will receive an error notification (no one lives forever)

If blood type is unknown, select A

Select "Detailed Detection Type"

Step 10b: Tap "Start Scan" button




Step 11b: Wait while MagnoSeek measures and analyzes.

Step 12b: Explore Your Results without Smart Phone



1



After the scan is finished, you will receive a notification that "detection is complete".

2


Your report overview will automatically generate.

3



Select "Reporting" at the top right corner to generate your overview report.

4




After the next window opens, select "collect".

5

Your detailed report will generate.

6



After the report is finished generating, select "Local".

Signal Interpretation Chart



Cells in a low-energy state
remain long-lived and stable



Optimal cellular regulation
is active and responsive



Cellular function aligns
with biological age



Mild dysfunction reduces
activity and vitality



Dysregulation significantly
decreases viability



Severe disruption leads to
functional decline and impairment



Warranty Coverage

Prife North America is dedicated to providing high-quality products. The MagnoSeek device comes with a two-year warranty from the date of purchase, while the MagnoSeek Wristband is covered for one year. This coverage offers peace of mind and reliable protection.

What This Warranty Covers:

- Manufacturing defects
- Workmanship issues
- Faulty materials

What This Warranty Does Not Cover:

- Damage due to misuse, accidents, or neglect
- Unauthorized repairs or modifications
- Normal wear and tear over time

How to Initiate a Warranty Claim:

Contact our Field Support team at support@prifemail.com or call 385-573-4931 with your order number and a description of the issue.

Ship the item to the address provided by our support team.

Warranty claims are typically processed within 10 business days after the item is received and inspected.

Return Shipping:

- Customers are responsible for return shipping costs unless the product is defective or the wrong item was received.
- We recommend using a trackable shipping service or purchasing shipping insurance for added peace of mind.

Return Shipping:

- Customers are responsible for return shipping costs unless the product is defective or the wrong item was received.
- We recommend using a trackable shipping service or purchasing shipping insurance for added peace of mind.



MagneSeek Support

MagneSeek User Manual.

Scan this code to access the official MagneSeek User Manual:



For any and all questions, comments, or concerns about MagneSeek, please don't hesitate to contact Prife North America Field Support.

Visit: PrifeShop.com to chat with our Field Support Team

Email: support@prife.com

Phone: (385) 573-4931

Serial Number#



Disclaimer: This system is not a medical device and does not diagnose, treat, or cure any condition. Results are for personal wellness reference only and should not replace professional medical advice, diagnosis, or treatment.